

SonoSite Selling Guide

M-Turbo 1.3

- Nerve exam setting now on the HFL50x
- Automated Cardiac Output from Doppler as licensed option on all phased array and the D2x
- L38xi with new housing and multi-angle biopsy brackets
- MSK exam type and worksheets on C60x, L25x, SLAx, HFL38x, L38xi and HFL50x
- Enhanced OB package with 42+ week measurements and expanded growth charts including Brenner table for EFW in the <10th percentile. Up to 10 GYN follicle measurements
- Vascular Volume flow units changed to ml/min and displayed on screen. Can angle correct on frozen image before saving. PW gate size can be used for diameter measurement if wanted